



## Nalgene® Pharma-Grade Platinum-Cured Silicone Tubing

### Product Images

---



### Description

---

Nalgene Pharma-Grade Platinum-Cured Silicone Tubing Ensure consistently reliable and accurate fluid transfers. These durable, high-purity tubings are ideal for pharmaceutical, laboratory and biotechnology applications. Certification/Compliance: Meets USP<sup>1</sup>87<sup>2</sup>gt; Cytotoxicity (Agar Diffusion and MEM Elution methods); USP<sup>1</sup>88<sup>2</sup>gt; Class VI; USP<sup>1</sup>66<sup>1</sup>gt; Physicochemical Tests; ISO10993-3 Bacterial Mutagenicity ndash; Ames Assay; ISO10993-4 Hemolysis Direct Contact, EP 3.1.9 Silicone Elastomer f Autoclavable: Autoclavable Inner Diameter (Anglo-American): 3/8 in. Outer Diameter (Anglo-American): 5/8 in. Overall length (anglo): 50 ft. Thickness (anglo): 1/8 in. Thickness (metric): 3.2 mm Type: Flexible Tubing Color: Translucent Hardness: Durometer 50(A) Flexibility: Bend radius 4x outside diameter Maximum pressure: 14 psig at 73 °F (23 °C) Unit Size: Case of 50

# Additional Information

|                      |         |
|----------------------|---------|
| No.                  | CD-2020 |
| Manufacturer (Brand) | Nalgene |

