



Nalgene® Pharma-Grade Platinum-Cured Silicone Tubing

Product Images



Description

Nalgene Pharma-Grade Platinum-Cured Silicone Tubing Ensure consistently reliable and accurate fluid transfers. These durable, high-purity tubings are ideal for pharmaceutical, laboratory and biotechnology applications. Certification/Compliance: Meets USP¹87²gt; Cytotoxicity (Agar Diffusion and MEM Elution methods); USP¹88²gt; Class VI; USP¹661²gt; Physicochemical Tests; ISO10993-3 Bacterial Mutagenicity ndash; Ames Assay; ISO10993-4 Hemolysis Direct Contact, EP 3.1.9 Silicone Elastomer f Autoclavable: Autoclavable Inner Diameter (Anglo-American): 1/4 in. Outer Diameter (Anglo): 1/2 in. Overall length (anglo): 50 ft. Thickness (Anglo): 1/8 in. Thickness (metric): 3.2 mm Type: Flexible Tubing Hardness: Durometer 50(A) Flexibility: Bend Radius 4xOD Maximum Pressure: 17 psig at 73deg;F (23deg;C) Unit Size: Case of 50

Additional Information

No.	CD-2016
Manufacturer (Brand)	Nalgene

