



Nalgene® Polycarbonate, End Cap

Product Images

Description

Nalgene Polycarbonate, End Cap Can be used with: Nalgene Balloon Bottle 2630 and 2261 Description: Two sanitary connectors with 0.75 in. (19 mm), polycarbonate inner diameter (anglo) tubing: 0.75 in. inner diameter (metric) tubing: 19 mm Type: End Cap Quantity per Case: 4 Quantity per Package: 1 For use with (devices): Nalgene Balloon Bottle 2630 and 2261 Unit Size: Case of 4

Additional Information

No.	CD-0696
Manufacturer (Brand)	Nalgene

