



Nalgene® Polycarbonate, Validation Bottle

Product Images

Description

Nalgene Polycarbonate, Validation Bottle Thermo Scientific Nalgene sample bottles are used for sample storage as part of quality assurance of high-value pharmaceuticals. Capacity (Anglo-American): 1 ounce Closure material: polypropylene Closure size: 20 - 415mm Inner diameter (metric) Neck: 14mm Type: Validation Bottle For use with (application): Perform validation testing Height (metric) with closure: 75mm Number per Case: 30 Number per Package: 30 Autoclavable Color: White (translucent) Mouth Style: Wide Unit Size: Case of 30

Additional Information

No.	CD-2112
Manufacturer (Brand)	Nalgene
VGKL number	228721271

