



Nalgene® Polycarbonate, Validation Bottle

Product Images

Description

Nalgene Polycarbonate, Validation Bottle Thermo Scientific Nalgene sample bottles are used for sample storage as part of quality assurance of high value pharmaceuticals. Capacity (Anglo): 64 ounces Closure Material: Polypropylene Closure Size: 53Bmm Inner Diameter (metric) Neck: 42mm Capacity (metric): 2L Type: Validation Bottle For use with (application): Perform validation testing Height (metric) with closure: 296mm Number per case: 12 Number per package: 12 Autoclavable Mouth Style: Wide Unit Size: Case of 12

Additional Information

No.	CD-2114
Manufacturer (Brand)	Nalgene

