



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene® polypropylene quick-connect taps

€396.00
plus VAT &
Shipping

Product Images



Description

Nalgene Polypropylene Quick Disconnect Cocks Can be used with: Nalgene Balloon Bottles Description: polypropylene quick disconnect cocks Inner diameter (anglo) tubing: 5/8 in. (Includes adapter for 1/4 or 5/16 in. tubing) Inner diameter (metric) tubing: 1.6 cm (Includes adapter for 0.6 or 0.8 cm tubing) Material: polypropylene, with HDPE tubing adapter Type: faucet Quantity per carton: 12 Quantity per package: 1 For use with (devices): Nalgene balloon bottles Unit Size: Case of 12

Additional Information

No.	CD-1773
Manufacturer (Brand)	Nalgene
VGKL number	228721623
Transport temperature	Room temperature

