



Nalgene® polypropylene vacuum flask

Product Images



Description

Nalgene™ Polypropylene Vacuum Flasks Replace glass flasks with Thermo Scientific™ Nalgene™ polypropylene vacuum flasks, molded in one piece for use with full vacuum and with angled tubing connection to prevent tip-over. Capacity (Anglo): 57.4 ounces Description: vacuum flask, polypropylene, 1700 ml Capacity (metric): 1700 ml Material: Polypropylene (PP) Quantity per carton: 1 Quantity per package: 1 Autoclavable: Yes Closure Size: No. 10 Stopper Unit Size: Each

Additional Information

No.	CD-2161
Manufacturer (Brand)	Nalgene

