



Nalgene® Polypropylene, End Cap

Product Images

Description

Nalgene Polypropylene, End Cap Can be used with: Nalgene Balloon Bottle 2630 and 2261 Description: Two sanitary 0.75 in. (19 mm), polypropylene inner diameter (Anglo) tubing: 0.75 in. inner diameter (metric) tubing: 19 mm Type: End Cap Quantity per Case: 4 Quantity per Package: 1 For use with (devices): Nalgene sanitary balloon bottles Unit Size: Case of 4

Additional Information

No.	CD-0697
Manufacturer (Brand)	Nalgene

