



Nalgene® Polypropylene, End Caps, mini

Product Images

Description

Nalgene Polypropylene, End Caps, mini Can be used with: Nalgene sanitary balloon bottles Cap Size: 19 mm Description: 0.75 in. Mini (19mm) Inner Diameter (metric) Tubing: 16 mm Type: Caps Quantity per carton: 1 Quantity per package: 1 For use with (devices): Nalgene Sanitary Balloon Bottles Unit Size: Each

Additional Information

No.	CD-0687
Manufacturer (Brand)	Nalgene

