



Nalgene® Rigid Containment Bottles

Product Images

Description

Nalgene Rigid Containment Bottles Containers developed specifically for the risk-averse Capacity: 5 L Capacity (metric): 5 L Capacity (anglo): 169 oz. Type: Bottle for Fluid Storage and Transport Lining Material: Silicone Certification/Compliance: USP Closure Liner: With Liner Closure Included: Yes Autoclavable: Not Autoclavable Closure Size: 48 mm Closure Material: Polypropylene Closure Type: Cap Graduated: Yes Height (Anglo-American): 11.77 in. Opening: Narrow Shape: Square Dispensing through closure: No Width (Anglo): 14.3 in. Unit Size: Case of 6

Additional Information

No.	CD-3652
Manufacturer (Brand)	Nalgene

