



Nalgene® Rigid Containment Bottles

Product Images

Description

Nalgene Rigid Containment Bottles Containers developed specifically for the risk-averse Capacity: 1 L Capacity (metric): 1 L Capacity (anglo): 33.8 oz. Type: Bottle for Fluid Storage and Transport Lining Material: Silicone Certification/Compliance: USP Class VI; USP 661; FDA 21CFR 177.1315 Closure Liner: With Liner Closure Included: Yes Autoclavable: Not Autoclavable Closure Size: 48 mm Closure Material: Polypropylene Closure Type: Cap Graduated: Yes Height (Anglo-American): 7.91 in. Opening: Narrow Shape: Square Dispensing through the closure: No Width (Anglo): 20.5 in. Unit Size: Case of,35

Additional Information

No.	CD-3649
Manufacturer (Brand)	Nalgene

