



Nalgene® tubes with linear barcode

Product Images



Description

Linear Barcode Tubes Simplify cryogenic tube tracking with linear barcodes printed directly on the side of the tubes. Single-handed

operation: Yes Capacity (Anglo-American): 0.0405 oz. Closure Material: Polyethylene Closure Type: Screw Cap DNase-free: Yes Package Quantity: 500/carton Graduated: No Type: Linear Barcode Certification/Compliance: Nalgene tubes are compliant with IATA for shipping. CE Marked for storage of diagnostic samples in the EU. Outer Diameter (Anglo-American): 0.5 inch Autoclavable: Not autoclavable Capacity (metric): 1.2 mL Thread Type: Outer Length (Anglo-American): 1.5 inches Number per carton: 500 Number per package: 25 RNase-free: yes Barcode: Yes Product Line: Nalgenetrade; Unit Size: Case of 500

Additional Information

No.	CD-1512
Manufacturer (Brand)	Nalgene

