



Nalgene® tubes with linear barcode

Product Images



Description

Linear Barcode Tubes Simplify cryogenic tube tracking with linear barcodes printed directly on the side of the tubes. Single-handed

operation: Yes Capacity (Anglo-American): 0.0676 oz. Closure Material: Polyethylene Closure Type: Screw cap DNase-free: Yes Package Quantity: 500/carton Graduated: No Type: Linear Barcode Certification/Compliance: Nunc tubes are IATA-approved for shipping. CE marked for storage of diagnostic specimens in EU countries. Description: Nalgene cryo vials with linear barcode (ink printer, code 128) Autoclavable: Non-autoclavable Capacity (metric): 2.0 ml Thread type: outer Number per carton: 500 Number per pack: 25 RNase-free: yes Barcode: Yes Product Line: Nunctrade; Unit Size: Case of 500

Additional Information

No.	CD-1513
Manufacturer (Brand)	Nalgene

