



Nalgene® Samco Tube Closures

Product Images



Description

Samco Tube Closures Depend on the reliable performance of these tube closures, available as flanged caps, two-position snap caps and universal caps. Description: Non-sterile universal cap for 12-16 mm tubes Quantity per box: 10000 Sterility: Non-sterile Type: Universal Plug Cap Unit Size: Case of 10,000 Description: Non-sterile universal cap for 12-16 mm tubes Quantity per box: 10000 Sterility: Non-sterile Type: Universal Plug Cap Unit Size: Case of 10,000 Description: Non-sterile universal cap for 12-16 mm tubes Quantity per box: 10000 Sterility: Non-sterile Type: Universal Plug Cap Unit Size: Case of 10,000

Additional Information

No.	CD-3666
Manufacturer (Brand)	Nalgene

