



Nalgene® sedimentation device kits

Product Images



Description

Nalgene Sedimentation Device Kits Thermo Scientific Nalgene Sedimentation Device Kits are ideal for the analysis of activated sludge samples. Certification/Compliance: Kit has been evaluated by the National Training and Operational Technology Center (NTOTC) of the U.S Closure Material: Polypropylene Description: ob Settlometer Number per box: 4 Number per package: 1 Type: Settlometer Height (Anglo): 7.5 in. Height (Anglo) Tubing: 4.7 in. Height (metric) Tube: 119mm Outer Diameter (Anglo): 5.2 in. Outer Diameter (Anglo) Tube: 0.6 in. Outer Diameter (Metric) Tube: 17mm Unit Size: Case of 4

Additional Information

No.	CD-0017
Manufacturer (Brand)	Nalgene
VGKL number	108710101

