



Nalgene® Silanized Disposable Culture Tubes

Product Images



Description

Silanized Disposable Culture Tubes Save time and minimize waste costs when performing quantitative analysis or storing materials with these ready-to-use, silanized disposable culture tubes. Description: Clear glass culture tubes Capacity (metric): 20 ml Material: Borosilicate glass (Type 1) Disposable: Yes Quantity: 250 pcs/pack, 1000 pcs/carton Type: Culture tubes Unit Size: Each

Additional Information

No.	CD-3663
Manufacturer (Brand)	Nalgene

