



## Nalgene® Silanized Disposable Culture Tubes

### Product Images

---



## Description

Silanized Disposable Culture Tubes Save time and minimize waste costs when performing quantitative analysis or storing materials with these ready-to-use, silanized disposable culture tubes. Description: Glass tube with screw thread Capacity (metric): 20 ml Closure: Screw cap (sold separately) Material: Borosilicate glass (Type 1) Disposable: Yes Quantity: 250 pcs/pack; 1000 pcs/carton Thread Size: 15-415 Type: Culture Tubes Unit Size: Each

## Additional Information

|                      |         |
|----------------------|---------|
| No.                  | CD-3717 |
| Manufacturer (Brand) | Nalgene |

