



Nalgene® Silanized Disposable Culture Tubes

Product Images



Description

Silanized Disposable Culture Tubes Save time and minimize waste costs when performing quantitative analysis or storing materials with these ready-to-use, silanized disposable culture tubes. Description: Screw thread glass tube Capacity (metric): 13 mL Closure: Screw Cap (Sold separately) Material: Borosilicate Glass (Type 1) Disposable: Yes Quantity: 250/pk, 1000/cs. Thread Size: 13-415 Type: Culture Tube Unit Size: Each

Additional Information

No.	CD-3716
Manufacturer (Brand)	Nalgene

