



Nalgene® Silanized Disposable Culture Tubes

Product Images



Description

Silanized Disposable Culture Tubes Save time and minimize waste costs when performing quantitative analysis or storing materials with these ready-to-use, silanized disposable culture tubes. Description: Clear glass culture tubes Capacity (metric): 20 mL Material: Borosilicate Glass (Type 1) Disposable Product: Yes Quantity: 250/pk, 1000/cs. Type: Culture Tube Unit Size: Each

Additional Information

No.	CD-3662
Manufacturer (Brand)	Nalgene

