



Nalgene® Silanized Glass Vials

Product Images



Description

Silanized Glass Vials Save time and minimize waste costs when performing quantitative analysis or for storing materials with these pre-treated, silanized vials. Description: amber glass vials 12 x 32 with threaded cap Outer diameter (anglo): 0.472 in. Color: Amber Volume (metric): 2 mL Height (anglo): 1.26 in. Lining Material: PTFE/Silicone Quantity Per Case: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Vial type: screw Cap type: screw cap with septum Quantity: 100 boxes Certification/Conformity: certificate of analysis (COA), certificate of compliance (COC) Type: vial Cap size: 8-425 Unit Size: Each Description: amber glass vial 12 x 32 with threaded cap Outer diameter (Anglo-American): 0.472 in. Color: Amber Volume (metric): 2 mL Height (anglo): 1.26 in. Lining Material: PTFE/Silicone Quantity Per Case: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Vial type: screw Cap type: screw cap with septum Quantity: 100 boxes Certification/Conformity: certificate of analysis (COA), certificate of compliance (COC) Type: vial Cap size: 8-425 Unit Size: Each Description: amber glass vial 12 x 32 with threaded cap Outer diameter (Anglo-American): 0.472 in. Color: Amber Volume (metric): 2 mL Height (anglo): 1.26 in. Lining Material: PTFE/Silicone Quantity Per Case: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Vial type: screw Cap type: screw cap with septum Quantity: 100 boxes Certification/Compliance: Certificate of Analysis (COA), Certificate of Conformity (COC) Type: vial Cap size: 8-425 Unit Size: Each

Additional Information

No.	CD-3694
Manufacturer (Brand)	Nalgene

