



Nalgene® Silanized Glass Vials

Product Images



Description

Silanized Glass Vials Save time and minimize waste costs when performing quantitative analysis or for storing materials with these pre-treated, silanized vials. Description: Amber glass vial with closure Outer diameter (anglo): 0.669 in. Volume (metric): 8 mL Height (anglo): 2.362 in. Lining material: PTFE/Foam Quantity per carton: 100 Volume (anglo): 0.27 oz. Total Volume (metric): 8 mL Vial Type: Screw Closure Type: Screw cap with septum Quantity: 100 boxes Certification/Conformance: Certificate of Analysis (COA), Certificate of Conformity (COC) Type: Vial Closure Size: 15-425 Unit Size: Each Description: Amber glass vial with closure Outer Diameter (Anglo-American): 0.669 in. Volume (metric): 8 mL Height (anglo): 2.362 in. Lining material: PTFE/Foam Quantity per carton: 100 Volume (anglo): 0.27 oz. Total Volume (metric): 8 mL Vial Type: Screw Closure Type: Screw cap with septum Quantity: 100 boxes Certification/Conformance: Certificate of Analysis (COA), Certificate of Conformity (COC) Type: Vial Closure Size: 15-425 Unit Size: Each Description: Amber glass vial with closure Outer Diameter (Anglo-American): 0.669 in. Volume (metric): 8 mL Height (anglo): 2.362 in. Lining material: PTFE/Foam Quantity per carton: 100 Volume (anglo): 0.27 oz. Total Volume (metric): 8 mL Vial Type: Screw Cap Closure Type: Screw Cap with Septum Quantity: 100 boxes Certification/Conformance: Certificate of Analysis (COA), Certificate of Conformity (COC) Type: Vial Cap Size: 15-425 Unit Size: Each

Additional Information

No.	CD-3665
Manufacturer (Brand)	Nalgene

