



Nalgene® Silanized Glass Vials

Product Images



Description

Silanized Glass Vials Save time and minimize waste costs when performing quantitative analysis or for storing materials with these pre-treated, silanized vials. Description: Clear glass vial 12 x 32 closed top closure Outer diameter (anglo): 0.472 in. Volume (metric): 2 mL Height (anglo): 1.26 in. Lining material: PTFE Quantity per box: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Closure material: PTFE Closure type: screw cap (closed top) Height (metric): 32 mm Quantity: 100 boxes Certification/Compliance: certificate of analysis (COA), certificate of compliance (COC) Type: vial Closure Size: 8-425 Unit Size: Each Description: clear glass vial 12 x 32 closed top closure Outer diameter (anglo): 0.472 in. Volume (metric): 2 mL Height (anglo): 1.26 in. Lining material: PTFE Quantity per box: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Closure material: PTFE Closure type: screw cap (closed top) Height (metric): 32 mm Quantity: 100 boxes Certification/Compliance: certificate of analysis (COA), certificate of compliance (COC) Type: vial Closure Size: 8-425 Unit Size: Each Description: clear glass vial 12 x 32 closed top closure Outer diameter (anglo): 0.472 in. Volume (metric): 2 mL Height (anglo): 1.26 in. Lining material: PTFE Quantity per box: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Closure material: PTFE Closure type: screw cap (closed top) Height (metric): 32 mm Quantity: 100 cartons Certification/Compliance: Certificate of Analysis (COA), Certificate of Conformity (COC) Type: Vial Closure Size: 8-425 Unit Size: Each

Additional Information

No.	CD-3698
Manufacturer (Brand)	Nalgene

