



Nalgene® Silanized Glass Vials

Product Images



Description

Silanized Glass Vials Save time and minimize waste costs when performing quantitative analysis or for storing materials with these pre-treated, silanized vials. Description: Amber glass vial 12 x 32 closed top closure Outer diameter (anglo): 0.472 in. Volume (metric): 2 mL Height (anglo): 1.26 in. Lining material: PTFE Quantity per box: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Vial Type: Screw Closure Material: PTFE Closure Type: Screw Cap (Closed-Top) Quantity: 100 Cases Certification/Conformance: Certificate of Analysis (COA), Certificate of Conformance (COC) Type: Vial Closure Size: 8-425 Unit Size: Each Description: Amber glass vial 12 x 32 closed top closure Outer diameter (Anglo-American): 0.472 in. Volume (metric): 2 mL Height (anglo): 1.26 in. Lining material: PTFE Quantity per box: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Vial Type: Screw Closure Material: PTFE Closure Type: Screw Cap (Closed-Top) Quantity: 100 Cases Certification/Conformance: Certificate of Analysis (COA), Certificate of Conformance (COC) Type: Vial Closure Size: 8-425 Unit Size: Each Description: Amber glass vial 12 x 32 closed top closure Outer diameter (Anglo-American): 0.472 in. Volume (metric): 2 mL Height (anglo): 1.26 in. Lining material: PTFE Quantity per box: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Vial Type: Screw Closure Material: PTFE Closure Type: Screw Cap (Closed-Top) Quantity: 100 Cases Certification/Compliance: Certificate of Analysis (COA), Certificate of Conformance (COC) Type: Vial Closure Size: 8-425 Unit Size: Each

Additional Information

No.	CD-3699
Manufacturer (Brand)	Nalgene

