



Nalgene® T-tube connectors

Product Images



Description

Nalgene T-Tube Connectors Maximize fluid flow with non-tapered T-type tube connectors that are made of polypropylene resin and are free of animal components. Description: T-type connector Autoclavable: Autoclavable Outer diameter (Anglo-American): 3/8 in. Overall length (anglo): 2.43 in. Overall length (metric): 61.72 mm Type: Tubing Connector For use with: Tubing I.D. 3/8 in. Package Quantity: 12 Unit Size: Case of 48

Additional Information

No.	CD-1712
Manufacturer (Brand)	Nalgene

