



Nalgene® table cooling container

Product Images



Description

Benchtop Cooling Container Protect enzymes, cells, reagents and solutions by maintaining freezing temperatures on the lab bench with the Benchtop Cooling Container. Array: 3 x 4 Construction/Design: gel-filled/ transparent lid Maintained temp, hours: 45 min Capacity (metric): 0.2 ml to 1.5 ml Holding capacity: 12 (cryotubes 0.2 to 1.5 ml) Temperature range: -20 °C Type: benchtop cooling container Height (metric): 125 mm Height (anglo): 4.94 in. Length (metric): 10.2 cm Length (anglo): 5.31 in. Width (metric): 14.4 cm Width (anglo): 4.25 in. Dimensions (L x W x H): 5-15/16 x 4-1/4 x 4-15/16 in. (151 x 108 x 125 mm) Description: For 12 0.2 - 1.5 ml

tubes, Tube arrangement: 3 x 4 Unit Size: Each

Additional Information

No.	CD-2167
Manufacturer (Brand)	Nalgene

