



Nalgene® Valved Barbed Quick Disconnects

Product Images



Description

Nalgene Valved Barbed Quick Disconnects Maintain the integrity of your fluids to be transferred by using Thermo Scientific Nalgene spout and valve quick disconnects. Description: Barbed, Valve Quick-Disconnect Inner Diameter (Anglo): 3/8 in. Overall length (anglo): 3.09 in. Overall length (metric): 78.5mm Type: Valved Barbed Quick Disconnect For use with: 3/8 in. (9.5mm) I.D. Package Quantity: 2/Pk., 12/Cs. Flexibility: Flexible Maximum Pressure: 2 psig Unit Size: Case of 12

Additional Information

No.	CD-1723
Manufacturer (Brand)	Nalgene

