



Nalgene® Valved Barbed Quick Disconnects

Product Images



Description

Nalgene Valved Barbed Quick Disconnects Maintain the integrity of your fluids to be transferred by using Thermo Scientific Nalgene Valved Barbed Quick Disconnects. Description: Barbed, Valve Quick-Disconnect Inside diameter (anglo): 0.25 in. Overall length (anglo): 3.28 in. Overall length (metric): 83.3 mm Type: Valved Barbed Quick Disconnect For use with: 1/4 in. (6.4mm) I.D. Packaging Quantity: 2/Pk., 12/Cs. Flexibility: Flexible Maximum Pressure: 2 psig Unit Size: Case of 12

Additional Information

No.	CD-1722
Manufacturer (Brand)	Nalgene

