



Nalgene® Two-Piece Buechner Funnel, Polypropylene

Product Images



16

Description

Two-Piece Nalgene Büchner Funnels, Polypropylene Thermo Scientific Nalgene Büchner Funnels made of polypropylene are lightweight and durable. Top diameter (metric): 93 mm Capacity (metric): 358 ml Autoclavable: Yes Description: Funnel Buchner PP 90 mm Length (metric) Stem: 90 mm Quantity: 4 Unit Size: Case of 4

Additional Information

No.	CD-1487
Manufacturer (Brand)	Nalgene

