



SI Analytics® Sterilin Autoclave Bag

Product Images

Description

Sterilin autoclave bags Sterilin autoclave bags are durable, disposable bags with a blue "biohazard" imprint. These bags are specifically designed for sterilization or disposal of contaminated waste in autoclaves or waste incinerators. Description: autoclave bags (HDPE) Material: HDPE (high density polyethylene) Quantity per carton: 200 Temperature range (metric): 121 deg;C Dimensions (L x W): 610 x 406 mm Length (anglo): 24 in. Width (anglo): 16 in. Unit Size: Case of 200

Additional Information

No.	SI-0231
Manufacturer (Brand)	SI Analytics

