



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## SI Analytics® Sterilin Autoclave Bag

**€350.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Sterilin autoclave bags Sterilin autoclave bags are durable, disposable bags with a blue "biohazard" imprint. These bags are specifically designed for sterilization or disposal of contaminated waste in autoclaves or waste incinerators. Description: autoclave bags (PP) Quantity per carton: 500- Temperature range (metric): 135 deg;C Dimensions (L x W): 601 x 406 mm Length (anglo): 24 in. Width (anglo): 16 in. Unit Size: Case of 500

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | SI-0236          |
| Manufacturer (Brand)  | SI Analytics     |
| Transport temperature | Room temperature |

