



SI Analytics® Sterilin Autoclave Bag

Product Images

Description

Sterilin autoclave bags Sterilin autoclave bags are durable, disposable bags with a blue "biohazard" imprint. These bags are specifically designed for sterilization or disposal of contaminated waste in autoclaves or waste incinerators. Description: autoclave bags (PP) Quantity per carton: 500 Temperature range (metric): 135deg;C Dimensions (L x W): 660 x 305 mm Length (anglo): 26 in. Width (anglo): 12 in. Unit Size: Case of 500

Additional Information

No.	SI-0234
Manufacturer (Brand)	SI Analytics

