



Thermo Scientific ABgene SmartCap Mats

Product Images



Description

ABgene SmartCap Mats The Thermo Scientific ABgene SmartCap Mat ensures a reliable seal from room temperatures as low as -80 C. Description: Seals 96 sample tubes simultaneously Type: sealing mat Number of wells: 96 Temperature range (metric): -80 °C For use with: Tubes 0.3 ml to 1.4 ml Unit Size: Each

Additional Information

No.	TA-0459
Manufacturer (Brand)	Thermo Scientific

