



Thermo Scientific Barnstead™ LabTower™ TII

Product Images

Description

Barnstead™ LabTower™ TII Convert tap water to high purity water with the Thermo Scientific™ Barnstead™ LabTower™ TII - a complete solution with treatment and tank. Frequency: 50/60 Hz Feedwater Connector: 3/4 in. NPT Type of Product Water: Type 2 TOC Monitor with UV Intensity Monitoring: No Type: Water Treatment System Electrical Resistance: 10 to 1 Mohm.cm Conductivity: 0.1 Rate 1.0µS/cm Feedwater Source: Cone Flow Rate: LabTower TII (UV) 40 with 3 l/min Pump Bacterial Content: 1 CFU/ml Particles 0.2 µm/ml: lt;1 Feedwater Monitoring: No Feedwater Pressure: 1 to 6 bar Operating Pressure: 2 to 6 bar Temperature Range (metric): 2 °C° to 35 °C° Bacterial retention: 99% Description: LabTower TII (UV) with 40 3L/min pump Permeate output: 40 l/h at 15 °C Electrical requirements: 100/240 V 50/60 Hz Unit Size: Each

Additional Information

No.	LD-1264
Manufacturer (Brand)	Thermo Scientific
VGKL number	108664640

