



Thermo Scientific CryoBox™ freezer racks

Product Images

Description

CryoBox™ Freezer Racks Conveniently store standard cryoboxes with these freezer racks made of durable stainless steel or aluminum that fit most standard freezers and chest freezers. Construction/Design: Vertical Dimensions (L x W x H): 5.625 x 26.5 x 5.531 in. (14.3 x 67 x 14 cm) Capacity: 12 CryoBoxes/Rack Type: Cryogenic Box Storage Rack, Vertical Unit Size: Each

Additional Information

No.	FE-0689
Manufacturer (Brand)	Thermo Scientific
VGKL number	290036700

