



# Thermo Scientific CryoMed Racks

## Product Images

## Description

CryoMed Racks For use with Thermo Scientific CryoMed controlled-rate freezers Description: 72 Packets ( 3 level) For use with (equipment): All CryoMed controlled-rate freezers Capacity: 72 Packets (3 level) Type: Skin -Rack Unit Size: Each

## Additional Information

|                      |                   |
|----------------------|-------------------|
| No.                  | LD-0762           |
| Manufacturer (Brand) | Thermo Scientific |

