



Thermo Scientific CryoMed Sensors

Product Images

Description

CryoMed Sensors For Thermo Scientific CryoMed freezers with step control. Description: Thermocouple Sensor Capacity (anglo): 7.06 cu. ft. Capacity (metric): 2 mL Connector type: male For use with (devices): For 1.2/2.0 cryo tube in CryoMed Type: Thermocouple Sensor Capacity: cryo tubes, jacketed straws, vials Unit Size: Each

Additional Information

No.	LD-0785
Manufacturer (Brand)	Thermo Scientific

