



Thermo Scientific CryoMed Sensors

Product Images

Description

CryoMed Sensors For Thermo Scientific CryoMed Step Control Freezers. Description: Thermocouple Sensor Capacity (anglo): 0.0001 cu. ft. Capacity (metric): 5 mL Connector Type: Male For use with (devices): For 4.0/5.0 cryo tube in CryoMed Type: Thermocouple Sensor Capacity: cryo tubes, jacketed straws, vials Unit Size: Each

Additional Information

No.	LD-0786
Manufacturer (Brand)	Thermo Scientific

