



Thermo Scientific F1-ClipTip™ GLP Kits

Product Images



Description

F1-ClipTip™ GLP Kits For maximum practicality, we have put together a range that includes a selection of F1-ClipTip pipettes of different volumes with compatible ClipTip pipette tips, plus a practical pipette stand for optimal performance of pipetting operations in the laboratory. Compatible tips: ClipTip Tips: 20, 200, 1000 Description: F1-ClipTip GLP Kit 1: 1-10micro;L, 10-100 micro;L and 100-1000 micro;L Number of Channels: 1 Volume (metric): 1 to 1000 micro;L Unit Size: Each

Additional Information

No.	PO-0247
Manufacturer (Brand)	Thermo Scientific

