



Thermo Scientific F1-ClipTip™ GLP Kits

Product Images

Description

F1-ClipTip™ GLP Kits For maximum practicality, we have put together a range that includes a selection of F1-ClipTip pipettes of different volumes with compatible ClipTip pipette tips, plus a practical pipette stand for optimal performance of pipetting operations in the laboratory. Compatible tips: ClipTip tips: 20, 200, 300, 1000 Description: F1-ClipTip GLP Kit 2: 2-20 micro;L, 20-200 micro;L, 100-1000 micro;L and 8-ch 30-300 micro;L Number of channels: 1 Volume (metric): 2 to 1000 and 30 to 300 micro;L Unit Size: Each

Additional Information

No.	PO-0249
Manufacturer (Brand)	Thermo Scientific

