



Thermo Scientific Adhesive Films

Product Images

Description

Adhesive Films Reduce your cost per plate in high-throughput applications with these adhesive films, which feature an automated format for efficient sealing. Description: Adhesive PCR Film - self-adhesive PCR film Material: polyethylene Type: Adhesive PCR Film - self-adhesive PCR film For use with (application): For sample storage, thermocycling with lid Adhesive: Yes Pierceable: Yes Temperature range (metric): -20 °F to +120 °C Plate Compatibility: 96-well plates, 384-well plates For use with: Seal-it 100 sealant Unit Size: Each

Additional Information

No.	TA-0469
Manufacturer (Brand)	Thermo Scientific

