



Thermo Scientific Cryoboxes

Product Images

Description

Cryoboxes For use with Thermo Scientific CryoPlus LN2 cryopreservation Description: 2.0 in., includes 40-cell divider Capacity (anglo): 0.0001 cu. ft. Capacity (metric): 5 mL For use with (devices): CryoPlus 1,2,3,4 Type: Fiberboard Box Capacity: Cryogenic Boxes Material: Stainless steel (box), polycarbonate (box), plastic Unit Size: Each

Additional Information

No.	LD-0741
Manufacturer (Brand)	Thermo Scientific

