



Thermo Scientific LN2 Hose

Product Images

Description

LN 2 Hose Ensure ease of use, safety and reliability in the event of a power failure or mechanical malfunction with this LN 2 hose.
Description: Excess hose (4 ft.) Temperature range: -150 °C to +1 °C For use with (devices): LNsub2/sub Supply Tank Defrost:
Automatic Type: LNsub2/sub Transfer Hose Unit Size: Each

Additional Information

No.	LD-0790
Manufacturer (Brand)	Thermo Scientific

