



Thermo Scientific Merlin™ Recirculating Chiller

Product Images

Description

Merlin™ Recirculating Chillers Save valuable space with Merlin's compact low-temperature chillers, which provide a continuous supply between 1000 and 4300 W of cooling capacity for critical applications. Description: M75 LR 230/60 T1 Pump Capacity: 3.5 @ 60 psi (13.0 @ 4.2 bar) Refrigerant: R134A Temperature Stability: plusmn;0.1deg;C Certification/Compliance: CE marked/UL listed Cooling Capacity: Up to 2,156 watts Electrical Requirements: 230 V 60 Hz Pump Type: T1 (Turbine Pump) Temperature Range: -15deg;C to +35deg;C Unit Size: Each

Additional Information

No.	LD-5448
Manufacturer (Brand)	Thermo Scientific

