



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Description: Creates a protective barrier against contamination Type: film, gasket or applicator Number of wells: 96 Number per box: 50 Number per package: 50 For use with (application): Aeration and storage of bacterial/cell cultures Product line: ABgene Adhesive: Gas permeable Sealant Puncturable: Not pierceable Temperature range (metric): -20 °C° to +80 °C° Plate Compatibility: All For use with: Used in conjunction with Thermo Scientific storage plates; 250050, 4110-11 Incubation Conditions: General Peelable: Yes Breathable: Yes

Additional Information

No.	AG-0012
Manufacturer (Brand)	Thermo Scientific

