



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Description: SEALING TAPE ADV POLYOLEFIN Material: polyethylene silicone Type: film, gasket or applicator Quantity per box: 100 Quantity per pack: 100 For use with (application): Q-PCR, microscopy, lowest fluorescence Product line: Nunc™ Adhesive: pressure sensitive Certified RNase/DNase-free: No Pierceable: No Temperature range (metric): -70 °C° to +100 °C° Plate Compatibility: 96-well plates, 384-well plates For use with: 250050, 4110-11 Solvent Or Alcohol Resistant: DMSO; ethanol; acetylnitrile; methanol; isopropanol Incubation Conditions: General Peelable: Yes Fluorescence: Low Autofluorescence Unit Size: Case of 100

Additional Information

No.	FE-0386
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900235307

