



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Description: Sealing Tape, Polyester, Clear, Roll Material: polyester/acrylic Type: film, seal or applicator Quantity per box: 1000 Quantity per package: 1000 For use with (application): Convenient roll dispenser, Multipurpose Product Line: Nunc™ Adhesive: Pressure Sensitive Certified RNase/DNase Free: Yes Pierceable: No Temperature range (metric): -20 °C° to +100 °C° Plate Compatibility: 96-well plates, 384-well plates Evaporation Barrier: Minimizes evaporation and cross contamination For use with: 250050, 4110-11 Incubation Conditions: General Peelable: Yes Unit Size: Case of 1000

Additional Information

No.	FE-0396
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900236707

