



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Description: White rayon sealing adhesive tape, RNase- and DNase-free Material: rayon/acrylate Type: film, seal or applicator Quantity per box: 50 Quantity per package: 20 For use with (application): Air permeable, Cell culture Product line: nunctrade; Adhesive: Pressure sensitive Certified RNase-/DNase-free: No Pierceable: Yes Temperature range (metric): -20 °C° to +80 °C° Plate Compatibility: 96-well plates, 384-well plates For use with: 250050, 4110-11 Incubation Conditions: General Peelable: Yes Breathable: Yes Unit Size: Case of 50

Additional Information

No.	FE-0406
Manufacturer (Brand)	Thermo Scientific

