



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Description: sealing film, acetate Material: acetate Type: film, seal or applicator Quantity per box: 100 For use with (application): For self-adhesive backing Product line: microtitertrade; Adhesive: Yes Temperature range (metric): -32 °C° to +93 °C° Plate Compatibility: All For use with: 250050, 4110-11 Incubation Conditions: General Peelable: Yes Unit Size: Each

Additional Information

No.	FE-0662
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900003501

