



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Color: Highly transparent Description: polyvinylidene chloride-coated polypropylene with acrylic adhesive Type: film, gasket or applicator Quantity per box: 100 Quantity per package: 100 For use with (application): For crystallography and genotyping Product line: ABgene Adhesive: Pressure sensitive Pierceable: Yes Temperature range (metric): -40 °C° to +80 °C° Plate Compatibility: All For use with: 250050, 4110-11 Incubation Conditions: General Peelable: Yes Unit Size: Case of 100

Additional Information

No.	AG-0007
Manufacturer (Brand)	Thermo Scientific

