



Thermo Scientific Nunc™ Sealing Tapes

Product Images

Description

Nunc Sealing Tapes Minimize evaporation and protect samples from contamination and leakage with Thermo Scientific Nunc Sealing Tape, ideal for all assay applications. Description: Sealing Tape Black Vinyl 50/cs Material: vinyl -acrylate Type: foil, gasket or applicator Quantity per box: 50 Quantity per package: 50 For use with (application): For lid or bottom of plate, photosensitive specimens Product line: Nunc™ Adhesive: Pressure sensitive Certified RNase/DNase free: No Pierceable: No Temperature range (metric): -40 °C° to +80 °C° Plate Compatibility: 96-well plates, 384-well plates For use with: 250050, 4110-11 Incubation Conditions: General Peelable: Yes Fluorescence: Low Autofluorescence Unit Size: Case of 50

Additional Information

No.	FE-0395
Manufacturer (Brand)	Thermo Scientific
VGKL number	2900236703

